4DryField® PH

PROVIDES HEMOSTASIS

**BLUTSTILLUNG
HÉMOSTASE
EMOSTASI
HEMOSTAZA
HEMOSTÁZA
HÆMOSTASE
HEMOSTAS
HEMOSTAZ**

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4DryField® PH – Hemostasis

Product Description

4DryField® PH is a medical device, an absorbable starch powder, for application to surgical wounds by physicians or medical-technical specialists.

4DryField® PH consists of hydrophilic microparticles produced from purified starch of plant-based origin. The free-flowing properties of the substance permit simple and safe application during surgical interventions. 4DryField® PH does not contain any components of animal or human origin. 4DryField® PH is sterilised, free of pyrogens, and biocompatible.

Mode of Action

Due to its molecular structure 4DryField® PH extracts water from the blood exuding from the wound. The cellular blood components (erythrocytes, platelets) and proteins (fibrinogen, thrombin) are concentrated, and therefore the coagulation cascade is accelerated.

4DryField® PH and blood form a hemostatic plug as a barrier against further blood leakage. Degradation of 4DryField® PH takes place via amylase and glucoamylase within days.

Indications

4DryField® PH is indicated in cases in which bleeding from arterioles, capillaries or veins during surgical interventions or in the case of injuries cannot or cannot effectively be controlled by conventional measures.

Contraindications

4DryField® PH is contraindicated in the presence of known intolerance to starch or starch-containing substances.

Caution: 4DryField® PH must not be administered intravascularly because this can lead to embolization and, in extreme cases, to death!

Administration of 4DryField® PH into the urinary system or into the respiratory tract is to be avoided because 4DryField® PH binds fluids and swells.

General Notes

The product must not be used if the package is damaged. Aseptic techniques should always be used. 4DryField® PH is only intended for use by physicians and medical-technical specialists. Therefore, these instructions for use cannot replace medical experience and judgement in the treatment of specific surgical situations.

Administration

1. Remove applicator cap with a slight bending and turning movement. 4DryField® PH is immediately ready for use.
2. Clean the bleeding area as thoroughly as possible (suction, standard gauze). 4DryField® PH displays its optimum effect if it is placed as closely to the source of bleeding as possible.
3. Administer 4DryField® PH liberally. The entire wound surface should be covered with 4DryField® PH. Wherever possible, contact of the applicator with blood or tissues is to be avoided in order to prevent blockage of the applicator.
4. Apply appropriate even pressure on the bleeding site, especially in the case of fairly large, parenchymatous bleeding, e.g. with standard gauze. Duration and magnitude of the pressure depend on the treated organ, severity of bleeding, and size of the wound surface. Before removing the standard gauze, if appropriate, rinse thoroughly with physiological saline. Bloody oozing may not require application of pressure.
5. If bloody oozing persists, remove excess 4DryField® PH by rinsing. Afterwards re-administer 4DryField® PH.
6. After contact with blood or other fluids 4DryField® PH swells within minutes to reach its maximum volume.
7. If 4DryField® PH is used in conjunction with an extracorporeal cardiopulmonary circulation (heart-lung machine) or autologous cell saving (blood processing) system, 40 µm filters and/or cell washing are to be employed.

Warnings

4DryField® PH is only intended for use by physicians or medical-technical specialists. 4DryField® PH does not replace thorough surgical hemostatic techniques via ligatures or other conventional measures. 4DryField® PH is not intended for primary treatment of coagulation disorders.

Excess 4DryField® PH is to be removed after hemostasis is reached in the area around the osseous foramina of the skull, in the area of the spinal column, of the spinal cord, or in the area of nerves in which nerve compression can occur (e.g. optic nerve, optic chiasm). In particular, 4DryField® PH that is still white should be removed. The likelihood of disturbances to normal functions or development of compression necroses brought on by swelling of 4DryField® PH is reduced upon removal of the excess.

Care should be taken to ensure that no compressions can be triggered when used in the area of the right atrium of the heart.

If applied in the cranial cavity one must take into account that 4DryField® PH swells due to the absorption of liquid. Introducing 4DryField® PH into cerebral fluid interspaces must be avoided.

Combination of 4DryField® PH with other hemostats or adhesives has not been investigated.

Polysaccharides such as 4DryField® PH can impair adhesion of acrylic-containing adhesives.

4DryField® PH has not been investigated in children or pregnant women. In new-borns up to ten months of age, amylase activity may be diminished, so that the absorption rate of products such as 4DryField® PH can be reduced.

4DryField® PH should not be used for control of postpartum bleeding or in cases of menorrhagia.

4DryField® PH should not be used if an infection is suspected. Signs of infection in an area in which 4DryField® PH was administered may require re-operation for insertion of a drainage.

So as to avoid contamination, do not reuse a unit of 4DryField® PH multiple times. Unused portions are to be disposed of.

Precautions

Particular diligence is required if 4DryField® PH is used in conjunction with a cardiopulmonary circulation or autologous cell saving system. Penetration of particles is prevented by use of 40 µm filters and cell washing. It is absolutely necessary that the instructions of the filter manufacturer are followed correctly. 4DryField® PH can only display its full effect in a dry state. Contact with saline or other solutions (e.g. antibiotic solutions) reduces the hemostatic potential.

4DryField® PH is supplied as a sterile product, cannot be re-sterilised because of loss of efficacy and, thus, cannot be re-used.

Other Information

In rare cases a selective increase in C-reactive protein (CRP) was reported during the first few days after the application of 4DryField® PH, which went down to levels within the normal range after a few days.

How Supplied

4DryField® PH can be supplied in bellows applicators containing 1 g, 3 g, 5 g, or 9 g. Applicators with 1 g are suitable for use in the dental field or ENT surgery. Applicators with 1 g, 3 g, 5 g, and 9 g are suitable for all other fields of surgery.

Storage

4DryField® PH should not be exposed to extreme temperatures. It is recommended to use 4DryField® PH up after opening the applicator. Unused portions are to be discarded in order to avoid contamination.

Disposal

The product is to be disposed according to the legal regulations for disposal of medical devices.

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Excess 4DryField® PH is to be removed when 4DryField® PH is used in the area around the osseous foramina of the skull, in the area of the spinal column, of the spinal cord, or in the area of nerves in which nerve compression can occur (e.g. optic nerve, optic chiasm). In particular, 4DryField® PH that is still white should be removed. The likelihood of disturbances to normal functions or development of compression necroses brought on by swelling of 4DryField® PH is reduced upon removal of the excess.

Care should be taken to ensure that no compressions can be triggered when used for adhesion prevention in the area of the right atrium of the heart.

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The safety and effectiveness of 4DryField® PH in combination with other medical devices regarding adhesion prevention has not been investigated.

Post-operative adhesions may occur even if 4DryField® PH is used. Possible causes may be insufficient hemostasis or improper use.

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Precautions

Particular diligence is required if 4DryField® PH is used in conjunction with a cardiopulmonary circulation or autologous cell saving system. Penetration of particles is prevented by use of 40 µm filters and cell washing. It is absolutely necessary that the instructions of the filter manufacturer are followed correctly.

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Mode of Action

The hydrophilic microparticles form a gel when fluid is absorbed due to the particular molecular structure of 4DryField® PH. For adhesion prevention, 4DryField® PH is transformed into a gel by sprinkling with isotonic saline solution. The gel functions as a temporary mechanical barrier that separates surgically traumatized mesothelial (or mesothelial-like) surfaces (areas that are at risk of the formation of adhesions). Degradation of 4DryField® PH as gel takes place via amylase and glucosylase within days.

Indications

4DryField® PH is indicated when the formation of post-operative adhesions is to be prevented following surgical interventions in cavities with a mesothelial (or mesothelial-like) lining.

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Administration

1. Remove applicator cap with a slight bending and turning movement. 4DryField® PH is immediately ready for use.
2. 4DryField® PH for adhesion prevention should only be employed in cavities with a mesothelial (or mesothelial-like) lining (e.g. abdominal cavity, chest cavity, pericardium). Care should be taken to ensure that there is no compression effect by 4DryField® PH as gel when body cavities are closed.
3. The affected areas should be cleaned as thoroughly as possible (suction, standard gauze) and rinsed.
4. Liberal application of 4DryField® PH. The entire defect should be liberally covered by 4DryField® PH. This also applies to large wound surfaces (i.e. pelvic exenteration, large peritoneal, pleural, pericardial defects). 4DryField® PH that is administered to non-traumatized mesothelial (or mesothelial-like) surfaces has no negative effects. Due to its outstanding biocompatibility, up to 1 g of 4DryField® PH per 1 kg body weight can be administered. Where possible, contact between the applicator and blood or tissue should be avoided, so as to prevent blockage of the applicator.
5. Following application of 4DryField® PH as powder, this is sprinkled with isotonic saline solution, until a gel completely covers the affected surface. White powder that has accidentally been left has no negative effects.
6. If 4DryField® PH is used in conjunction with an extracorporeal cardiopulmonary circulation (heart-lung machine) or autologous cell saving (blood processing) system, 40 µm filters and/or cell washing are to be employed.

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© PlantTec Medical GmbH 2015
4DryField® PH and design is a trademark of
PlantTec Medical GmbH, Bad Bevensen, Germany